




# Bio-bin® Fold Flat Instruction sheet



1. Take flat Bio-bin® out of box and hold as shown




2. Pull Bio-bin® into shape



3. Fold sides outwards to prevent bin folding flat again




4. Bio-bin® is ready to use for non-sharp infectious waste




5. Fold top down towards centre



6. Press firmly in the centre of flap to temporary seal



7. For a permanent seal, peel off release paper



8. Press together firmly and then fold down



9. Stick tabs to sides of Bio-bin®, lining up with printed area.